## High Efficiency Direct Connect (DC) Spray Chamber



- · No O-rings.
- DC connection is inert all PEEK construction.
- Eliminates ball joint clamps that rust over time, causing lab contamination.
- Guarantees perfect alignment each and every time for improved precision and better transport efficiency.
- Tracey<sup>™</sup> DC Spray chamber design is available in Glass, PEEK, PFA and Quartz material.
- Low volume (30mL) cyclonic spray chamber for faster washout with Helix CT technology.

Scan the QR Code to watch the DC Spray Chamber Installation video:



